


TAKE YOUR OBAGI[®] SKINCARE RESULTS TO THE NEXT LEVEL WITH A SKINTRINSIQ[™] TREATMENT

EXTRACT




Built-up impurities are removed from the skin's surface.

INFUSE




Proven Obagi products are rapidly distributed to absorb into the skin.

LOCK



Once infusion is complete, product remains in the skin to continue working after the treatment.

LED BENEFITS



The Skintrinsiq device has the unique ability to integrate both blue and red LED light without additional treatment time.

